

# ABHD12B Polyclonal Antibody

Cat No: HR1AP5971

For research use only

## Overview

Product Name	ABHD12B Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	ABHD12B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABHD12B; C14orf29; Abhydrolase domain-containing protein 12B
Background	similarity:Belongs to the serine esterase family.,