

ACSS1 Polyclonal Antibody

Cat No: HR1AP9639

For research use only

Overview

Product Name	ACSS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	ACSS1 ACAS2L KIAA1846
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACSS1 ACAS2L KIAA1846
Background	acyl-CoA synthetase short-chain family member 1(ACSS1) Homo sapiens This gene encodes a mitochondrial acetyl-CoA synthetase enzyme. A similar protein in mice plays an important role in the tricarboxylic acid cycle by catalyzing the conversion of acetate to acetyl CoA. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011],