

UBA1 Polyclonal Antibody

Cat No: HR1AP5334

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | UBA1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 118kDa |
| GeneID?Human? | UBA1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | UBA1; A1S9T; UBE1; Ubiquitin-like modifier-activating enzyme 1; Protein A1S9; Ubiquitin-activating enzyme E1 |
| Background | ubiquitin like modifier activating enzyme 1(UBA1) Homo sapiens The protein encoded by this gene catalyzes the first step in ubiquitin conjugation to mark cellular proteins for degradation. This gene complements an X-linked mouse temperature-sensitive defect in DNA synthesis, and thus may function in DNA repair. It is part of a gene cluster on chromosome Xp11.23. Alternatively spliced transcript variants that encode the same protein have been described. [provided by RefSeq, Jul 2008]. |