

ZADH2 Polyclonal Antibody

Cat No: HR1AP4833

For research use only

Overview

Product Name	ZADH2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	ZADH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZADH2; Zinc-binding alcohol dehydrogenase domain-containing protein 2
Background	similarity:Belongs to the zinc-containing alcohol dehydrogenase family. Quinone oxidoreductase subfamily.,