

UBE2T Polyclonal Antibody

Cat No: HR1AP4777

For research use only

Overview

Product Name	UBE2T Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	UBE2T
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	UBE2T; HSPC150; PIG50; Ubiquitin-conjugating enzyme E2 T; Cell proliferation-inducing gene 50 protein; Ubiquitin carrier protein T; Ubiquitin-protein ligase T
Background	ubiquitin conjugating enzyme E2 T(UBE2T) Homo sapiens The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2015],