

PFAAP5 Polyclonal Antibody

Cat No: HR1AP5583

For research use only

Overview

Product Name	PFAAP5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	N4BP2L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N4BP2L2; CG005; PFAAP5; NEDD4-binding protein 2-like 2; Phosphonoformate immuno-associated protein 5
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,