

PC-PLD4 Polyclonal Antibody

Cat No: HR1AP5812

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | PC-PLD4 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 50kDa |
| GeneID?Human? | PLD4 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | PLD4; C14orf175; Phospholipase D4; PLD 4; Choline phosphatase 4; Phosphatidylcholine-hydrolyzing phospholipase D4 |
| Background | catalytic activity:A phosphatidylcholine + H(2)O = choline + a phosphatidate.,similarity:Belongs to the phospholipase D family.,similarity:Contains 2 PLD phosphodiesterase domains., |