

PAH Polyclonal Antibody

Cat No: HR1AP4253

For research use only

Overview

Product Name	PAH Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	PAH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PAH; Phenylalanine-4-hydroxylase; PAH; Phe-4-monooxygenase
Background	phenylalanine hydroxylase(PAH) Homo sapiens PAH encodes the enzyme phenylalanine hydroxylase that is the rate-limiting step in phenylalanine catabolism. Deficiency of this enzyme activity results in the autosomal recessive disorder phenylketonuria. [provided by RefSeq, Jul 2008],