

# Neuralized-2 Polyclonal Antibody

Cat No: HR1AP4036

For research use only

## Overview

Product Name	Neuralized-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	NEURL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NEURL2; C20orf163; Neuralized-like protein 2
Background	neuralized E3 ubiquitin protein ligase 2(NEURL2) Homo sapiens This gene encodes a protein that is involved in the regulation of myofibril organization. This protein is likely the adaptor component of the E3 ubiquitin ligase complex in striated muscle, and it regulates the ubiquitin-mediated degradation of beta-catenin during myogenesis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2013],