

GUSBP1 Polyclonal Antibody

Cat No: HR1AP6026

For research use only

Overview

Product Name	GUSBP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	GUSBP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GUSBP1; SMA3; Putative beta-glucuronidase-like protein SMA3
Background	Binds GTP, GDP and GMP.