

GM2/GD2 synthase Polyclonal Antibody

Cat No: HR1AP3551

For research use only

Overview

Product Name	GM2/GD2 synthase Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	B4GALNT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	B4GALNT1; GALGT; SIAT2; Beta-1; 4 N-acetylgalactosaminyltransferase 1; (N-acetylneuraminyl)-galactosylglucosylceramide; GM2/GD2 synthase; GalNAc-T
Background	beta-1,4-N-acetyl-galactosaminyltransferase 1(B4GALNT1) Homo sapiens GM2 and GD2 gangliosides are sialic acid-containing glycosphingolipids. GalNAc-T is the enzyme involved in the biosynthesis of G(M2) and G(D2) glycosphingolipids. GalNAc-T catalyzes the transfer of GalNAc into G(M3) and G(D3) by a beta-1,4 linkage, resulting in the synthesis of G(M2) and G(D2), respectively. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2013],