

CYP2W1 Polyclonal Antibody

Cat No: HR1AP3218

For research use only

Overview

Product Name	CYP2W1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	CYP2W1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYP2W1; Cytochrome P450 2W1; CYP11W1
Background	cytochrome P450 family 2 subfamily W member 1(CYP2W1) Homo sapiens This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. [provided by RefSeq, Jul 2008],