

CYP2U1 Polyclonal Antibody

Cat No: HR1AP3217

For research use only

Overview

Product Name	CYP2U1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	CYP2U1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYP2U1; Cytochrome P450 2U1
Background	cytochrome P450 family 2 subfamily U member 1(CYP2U1) Homo sapiens This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This enzyme is a hydroxylase that metabolizes arachidonic acid, docosahexaenoic acid, and other long chain fatty acids. [provided by RefSeq, Jul 2008],