

CA XIII Polyclonal Antibody

Cat No: HR1AP2926

For research use only

Overview

Product Name	CA XIII Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	CA13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CA13; Carbonic anhydrase 13; Carbonate dehydratase XIII; Carbonic anhydrase XIII; CA-XIII
Background	carbonic anhydrase 13(CA13) Homo sapiens Carbonic anhydrases (CAs) are a family of zinc metalloenzymes. For background information on the CA family, see MIM 114800.[supplied by OMIM, Mar 2008],