

FAT10 Polyclonal Antibody

Cat No: HR1AP9471

For research use only

Overview

Product Name	FAT10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	UBD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	UBD; FAT10; Ubiquitin D; Diubiquitin; Ubiquitin-like protein FAT10
Background	similarity:Contains 2 ubiquitin-like domains.,subunit:Interacts with MAD2.,