

N4BP1 Polyclonal Antibody

Cat No: HR1AP3999

For research use only

Overview

Product Name	N4BP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	N4BP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N4BP1; KIAA0615; NEDD4-binding protein 1; N4BP1
Background	PTM:Monoubiquitinated.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the N4BP1 family.,subcellular location:In discrete circular structures.,subunit:Binds NEDD4.,tissue specificity:Detected in heart, lung, brain, liver, skeletal muscle, pancreas, kidney, spleen, testis and ovary.,