

## Calpain 12 Polyclonal Antibody

Cat No: HR1AP2943

For research use only

### Overview

Product Name	Calpain 12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	CAPN12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CAPN12; Calpain-12; Calcium-activated neutral proteinase 12; CANP 12
Background	calpain 12(CAPN12) Homo sapiens The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes a member of the calpain large subunit family. [provided by RefSeq, Jun 2012].