

UBE4A Polyclonal Antibody

Cat No: HR1AP4921

For research use only

Overview

Product Name	UBE4A Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	122kDa
GeneID?Human?	UBE4A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	UBE4A; KIAA0126; Ubiquitin conjugation factor E4 A
Background	<p>ubiquitination factor E4A(UBE4A) Homo sapiens This gene encodes a member of the U-box ubiquitin ligase family. The encoded protein is involved in multiubiquitin chain assembly and plays a critical role in chromosome condensation and separation through the polyubiquitination of securin. Autoantibodies against the encoded protein may be markers for scleroderma and Crohn's disease. A pseudogene of this gene is located on the long arm of chromosome 3. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Feb 2011],</p>