

CYP4Z1 Polyclonal Antibody

Cat No: HR1AP6201

For research use only

Overview

Product Name	CYP4Z1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CYP4Z1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYP4Z1; Cytochrome P450 4Z1; CYP4Z1
Background	CYP4Z1 encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. CYP4Z1 is part of a cluster of cytochrome P450 genes on chromosome 1p33.