

ACOT9 Polyclonal Antibody

Cat No: HR1AP6480

For research use only

Overview

Product Name	ACOT9 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	ACOT9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACOT9; CGI-16; Acyl-coenzyme A thioesterase 9; mitochondrial; Acyl-CoA thioesterase 9; Acyl-CoA thioester hydrolase 9
Background	acyl-CoA thioesterase 9(ACOT9) Homo sapiens The protein encoded by this gene is a mitochondrial acyl-CoA thioesterase of unknown function. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010],