

# Hepatic Lipase Polyclonal Antibody

Cat No: HR1AP7230

For research use only

## Overview

Product Name	Hepatic Lipase Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	LIPC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LIPC; HTGL; Hepatic triacylglycerol lipase; HL; Hepatic lipase; Lipase member C
Background	lipase C, hepatic type(LIPC) Homo sapiens LIPC encodes hepatic triglyceride lipase, which is expressed in liver. LIPC has the dual functions of triglyceride hydrolase and ligand/bridging factor for receptor-mediated lipoprotein uptake. [provided by RefSeq, Jul 2008],