

HCCS Polyclonal Antibody

Cat No: HR1AP3607

For research use only

Overview

Product Name	HCCS Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	HCCS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HCCS; CCHL; Cytochrome c-type heme lyase; CCHL; Holocytochrome c-type synthase
Background	holocytochrome c synthase(HCCS) Homo sapiens The protein encoded by this gene is an enzyme that covalently links a heme group to the apoprotein of cytochrome c. Defects in this gene are a cause of microphthalmia syndromic type 7 (MCOPS7). Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jan 2010],