

## DHODH Polyclonal Antibody

Cat No: HR1AP6116

For research use only

### Overview

Product Name	DHODH Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	DHODH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DHODH; Dihydroorotate dehydrogenase; quinone); mitochondrial; DHOdehase; Dihydroorotate oxidase
Background	dihydroorotate dehydrogenase (quinone)(DHODH) Homo sapiens The protein encoded by this gene catalyzes the fourth enzymatic step, the ubiquinone-mediated oxidation of dihydroorotate to orotate, in de novo pyrimidine biosynthesis. This protein is a mitochondrial protein located on the outer surface of the inner mitochondrial membrane. [provided by RefSeq, Jul 2008],