

# Cathepsin E Polyclonal Antibody

Cat No: HR1AP2982

For research use only

## Overview

Product Name	Cathepsin E Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	CTSE
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CTSE; Cathepsin E
Background	<p>cathepsin E(CTSE) Homo sapiens This gene encodes a member of the A1 family of peptidases. Alternative splicing of this gene results in multiple transcript variants. At least one of these variants encodes a preproprotein that is proteolytically processed to generate the mature enzyme. This enzyme, an aspartic endopeptidase, may be involved in antigen processing and the maturation of secretory proteins. Elevated expression of this gene has been observed in neurodegeneration. [provided by RefSeq, Nov 2015],</p>