

AMPD2 Polyclonal Antibody

Cat No: HR1AP2757

For research use only

Overview

Product Name	AMPD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	AMPD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AMPD2; AMP deaminase 2; AMP deaminase isoform L
Background	adenosine monophosphate deaminase 2(AMPD2) Homo sapiens The protein encoded by this gene is important in purine metabolism by converting AMP to IMP. The encoded protein, which acts as a homotetramer, is one of three AMP deaminases found in mammals. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012],