

NDUFA4L2 Polyclonal Antibody

Cat No: HR1AP4018

For research use only

Overview

Product Name	NDUFA4L2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	NDUFA4L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NDUFA4L2; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2; NADH-ubiquinone oxidoreductase MLRQ subunit homolog; NUOMS
Background	similarity:Belongs to the complex I NDUFA4 subunit family.,