

Atg16 Polyclonal Antibody

Cat No: HR1AP7811

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | Atg16 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 70kDa |
| GeneID?Human? | ATG16L1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | ATG16L1; APG16L; Autophagy-related protein 16-1; APG16-like 1 |
| Background | <p>autophagy related 16 like 1(ATG16L1) Homo sapiens The protein encoded by this gene is part of a large protein complex that is necessary for autophagy, the major process by which intracellular components are targeted to lysosomes for degradation. Defects in this gene are a cause of susceptibility to inflammatory bowel disease type 10 (IBD10). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jun 2010],</p> |