

# ABHD11 Polyclonal Antibody

Cat No: HR1AP2681

For research use only

## Overview

Product Name	ABHD11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	ABHD11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABHD11; WBSR21; PP1226; Abhydrolase domain-containing protein 11; Williams-Beuren syndrome chromosomal region 21 protein
Background	abhydrolase domain containing 11(ABHD11) Homo sapiens This gene encodes a protein containing an alpha/beta hydrolase fold domain. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. [provided by RefSeq, Mar 2016],