

# FAM35A/B Polyclonal Antibody

Cat No: HR1AP7773

For research use only

## Overview

Product Name	FAM35A/B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	FAM35A/FAM35B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FAM35A; Protein FAM35A; FAM35B; Protein FAM35B
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the FAM35 family.,