

SENP1 Polyclonal Antibody

Cat No: HR1AP4530

For research use only

Overview

Product Name	SENP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	SENP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SENP1; Sentrin-specific protease 1; Sentrin/SUMO-specific protease SENP1
Background	SUMO1/sentrin specific peptidase 1(SENP1) Homo sapiens This gene encodes a cysteine protease that specifically targets members of the small ubiquitin-like modifier (SUMO) protein family. This protease regulates SUMO pathways by deconjugating sumoylated proteins. This protease also functions to process the precursor SUMO proteins into their mature form. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012],