

IFITM2 Polyclonal Antibody

Cat No: HR1AP5373

For research use only

Overview

Product Name	IFITM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	IFITM2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IFITM2; Interferon-induced transmembrane protein 2; Dispanin subfamily A member 2c; DSPA2c; Interferon-inducible protein 1-8D
Background	induction:By alpha and gamma interferons.,similarity:Belongs to the CD225 family.,