

NDUFS3 Polyclonal Antibody

Cat No: HR1AP7461

For research use only

Overview

Product Name	NDUFS3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	NDUFS3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NDUFS3; NADH dehydrogenase [ubiquinone] iron-sulfur protein 3; mitochondrial; Complex I-30kD; CI-30kD; NADH-ubiquinone oxidoreductase 30 kDa subunit
Background	NADH:ubiquinone oxidoreductase core subunit S3(NDUFS3) Homo sapiens This gene encodes one of the iron-sulfur protein (IP) components of mitochondrial NADH:ubiquinone oxidoreductase (complex I). Mutations in this gene are associated with Leigh syndrome resulting from mitochondrial complex I deficiency.[provided by RefSeq, Apr 2009],