

# ACC? (phospho Ser80) Polyclonal Antibody

Cat No: HR1AP2578

For research use only

## Overview

Product Name	ACC? (phospho Ser80) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	280kDa
GeneID?Human?	ACACA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACACA; ACAC; ACC1; ACCA; Acetyl-CoA carboxylase 1; ACC1; ACC-alpha
Background	<p>acetyl-CoA carboxylase alpha(ACACA) Homo sapiens Acetyl-CoA carboxylase (ACC) is a complex multifunctional enzyme system. ACC is a biotin-containing enzyme which catalyzes the carboxylation of acetyl-CoA to malonyl-CoA, the rate-limiting step in fatty acid synthesis. There are two ACC forms, alpha and beta, encoded by two different genes. ACC-alpha is highly enriched in lipogenic tissues. The enzyme is under long term control at the transcriptional and translational levels and under short term regulation by the phosphorylation/dephosphorylation of targeted serine residues and by allosteric transformation by citrate or palmitoyl-CoA. Multiple alternatively spliced transcript variants divergent in the sequence and encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>