



USP30 Polyclonal Antibody

Cat No: HR1AP4789

For research use only

Overview

Product Name	USP30 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kd,unspecial bands(65kd,105kd,130kd)in Hela and HCT116 cell.kDa
GeneID?Human?	USP30
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	USP30; Ubiquitin carboxyl-terminal hydrolase 30; Deubiquitinating enzyme 30; Ubiquitin thioesterase 30; Ubiquitin-specific-processing protease 30; Ub-specific protease 30
Background	ubiquitin specific peptidase 30(USP30) Homo sapiens USP30, a member of the ubiquitin-specific protease family (see USP1, MIM 603478), is a novel mitochondrial deubiquitinating (DUB) enzyme (Nakamura and Hirose, 2008 [PubMed 18287522]).[supplied by OMIM, Dec 2008],