

Mox1 Polyclonal Antibody

Cat No: HR1AP9746

For research use only

Overview

Product Name	Mox1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	NOX1 MOX1 NOH1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NADPH oxidase 1 (NOX-1) (EC 1.-.-) (Mitogenic oxidase 1) (MOX-1) (NADH/NADPH mitogenic oxidase subunit P65-MOX) (NOH-1)
Background	NADPH oxidase 1(NOX1) Homo sapiens This gene encodes a member of the NADPH oxidase family of enzymes responsible for the catalytic one-electron transfer of oxygen to generate superoxide or hydrogen peroxide. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2012],