

SPZ1 Polyclonal Antibody

Cat No: HR1AP8882

For research use only

Overview

Product Name	SPZ1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	133kDa
GeneID?Human?	SPZ1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SPZ1; TSP1; Spermatogenic leucine zipper protein 1; Testis-specific protein 1; Testis-specific protein NYD-TSP1
Background	spermatogenic leucine zipper 1(SPZ1) Homo sapiens This gene encodes a bHLH-zip transcription factor which functions in the mitogen-activate protein kinase (MAPK) signaling pathway. Because of its role in the upregulation of cell proliferation and tumorigenesis, this gene may serve as a target for Ras-induced tumor treatments. [provided by RefSeq, Oct 2011],