

## S-100A16 Polyclonal Antibody

Cat No: HR1AP5933

For research use only

### Overview

Product Name	S-100A16 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	S100A16
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	S100A16; S100F; AAG13; Protein S100-A16; Aging-associated gene 13 protein; Protein S100-F; S100 calcium-binding protein A16
Background	similarity:Belongs to the S-100 family.,similarity:Contains 2 EF-hand domains.,