

# POTE-8 Polyclonal Antibody

Cat No: HR1AP6928

For research use only

## Overview

Product Name	POTE-8 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	POTEA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	POTEA; A26A1; POTE8; POTE ankyrin domain family member A; ANKRD26-like family A member 1; Prostate; ovary, testis-expressed protein on chromosome 8; POTE-8
Background	similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,