

Dynein LC 4 Polyclonal Antibody

Cat No: HR1AP5537

For research use only

Overview

Product Name	Dynein LC 4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	DNAL4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DNAL4; Dynein light chain 4; axonemal
Background	dynein axonemal light chain 4(DNAL4) Homo sapiens This gene encodes an axonemal dynein light chain which functions as a component of the outer dynein arms complex. This complex acts as the molecular motor that provides the force to move cilia in an ATP-dependent manner. The encoded protein is expressed in tissues with motile cilia or flagella and may be involved in the movement of sperm flagella. [provided by RefSeq, Dec 2014],