

## Dynein IC2 Polyclonal Antibody

Cat No: HR1AP8288

For research use only

### Overview

Product Name	Dynein IC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	DNAI2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DNAI2; Dynein intermediate chain 2; axonemal; Axonemal dynein intermediate chain 2
Background	dynein axonemal intermediate chain 2(DNAI2) Homo sapiens The protein encoded by this gene belongs to the dynein intermediate chain family, and is part of the dynein complex of respiratory cilia and sperm flagella. Mutations in this gene are associated with primary ciliary dyskinesia type 9. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Mar 2010],