

Ameloblastin Polyclonal Antibody

Cat No: HR1AP6550

For research use only

Overview

Product Name	Ameloblastin Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	AMBN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AMBN; Ameloblastin
Background	ameloblastin(AMBN) Homo sapiens This gene encodes the nonamelogenin enamel matrix protein ameloblastin. The encoded protein may be important in enamel matrix formation and mineralization. This gene is located in the calcium-binding phosphoprotein gene cluster on chromosome 4. Mutations in this gene may be associated with dentinogenesis imperfect and autosomal dominant amylogenesis imperfect. [provided by RefSeq, Aug 2011],