

S-100Z Polyclonal Antibody

Cat No: HR1AP6110

For research use only

Overview

Product Name	S-100Z Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	S100Z
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	S100Z; Protein S100-Z; S100 calcium-binding protein Z
Background	S100 calcium binding protein Z(S100Z) Homo sapiens Members of the S100 protein family contain 2 calcium-binding EF-hands and exhibit cell-type specific expression patterns. For additional background information on S100 proteins, see MIM 114085.[supplied by OMIM, Mar 2008],