

# POTE14/22 Polyclonal Antibody

Cat No: HR1AP4357

For research use only

## Overview

|                       |   |
|-----------------------|---|
| Product Name          | POTE14/22 Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,ELISA  |
| Species Reactivity    | Human   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 67kDa   |
| GeneID?Human?         | POTEG   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | POTEG; A26C2; POTE14; POTE ankyrin domain family member G; ANKRD26-like family C member 2; Prostate; ovary, testis-expressed protein on chromosome 14; POTE-14; POTEH; A26C3; ACTBL1; POTE22; POTE anky |
| Background            | similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,  |