

## CLIP4 Polyclonal Antibody

Cat No: HR1AP3111

For research use only

### Overview

Product Name	CLIP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	CLIP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLIP4; RSNL2; CAP-Gly domain-containing linker protein 4; Restin-like protein 2
Background	similarity:Contains 3 ANK repeats.,similarity:Contains 3 CAP-Gly domains.,