

Centrobin Polyclonal Antibody

Cat No: HR1AP3059

For research use only

Overview

Product Name	Centrobin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	CNTROB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNTROB; LIP8; PP1221; Centrobin; Centrosomal BRCA2-interacting protein; LYST-interacting protein 8
Background	centrobin, centriole duplication and spindle assembly protein(CNTROB) Homo sapiens This gene encodes a centrosomal protein that interacts with BRCA2, and is required for centriole duplication and cytokinesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011],