

## CaBP7 Polyclonal Antibody

Cat No: HR1AP6091

For research use only

### Overview

Product Name	CaBP7 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	CABP7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CABP7; CALN2; Calcium-binding protein 7; CaBP7; Calneuron II; Calneuron-2
Background	similarity:Contains 2 EF-hand domains.,