

## CLIP-170 Polyclonal Antibody

Cat No: HR1AP3109

For research use only

### Overview

Product Name	CLIP-170 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	161kDa
GeneID?Human?	CLIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLIP1; CYLN1; RSN; CAP-Gly domain-containing linker protein 1; Cytoplasmic linker protein 1; Cytoplasmic linker protein 170 alpha-2; CLIP-170; Reed-Sternberg intermediate filament-associated protein;
Background	CAP-Gly domain containing linker protein 1(CLIP1) Homo sapiens The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],