

Arp3 Polyclonal Antibody

Cat No: HR1AP2818

For research use only

Overview

Product Name	Arp3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	ACTR3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACTR3; ARP3; Actin-related protein 3; Actin-like protein 3
Background	<p>ARP3 actin related protein 3 homolog(ACTR3) Homo sapiens The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Mar 2013],</p>