

CABP2 Polyclonal Antibody

Cat No: HR1AP9406

For research use only

Overview

Product Name	CABP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	CABP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Calcium-binding protein 2 (CaBP2)
Background	<p>calcium binding protein 2(CABP2) Homo sapiens This gene belongs to a subfamily of calcium binding proteins that share similarity to calmodulin. Like calmodulin, these family members can likely stimulate calmodulin-dependent kinase II and the protein phosphatase calcineurin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016],</p>