

## Desmin (phospho Ser60) Polyclonal Antibody

Cat No: HR1AP6103

For research use only

### Overview

Product Name	Desmin (phospho Ser60) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	DES
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DES; Desmin
Background	desmin(DES) Homo sapiens This gene encodes a muscle-specific class III intermediate filament. Homopolymers of this protein form a stable intracytoplasmic filamentous network connecting myofibrils to each other and to the plasma membrane. Mutations in this gene are associated with desmin-related myopathy, a familial cardiac and skeletal myopathy (CSM), and with distal myopathies. [provided by RefSeq, Jul 2008],